




SU HANDLE

SY HANDLE

RD HANDLE

RY HANDLE



Catheter Diameter

Catheter Length

Basket Opening



STERILE EO

Coding

BKH 5 20 12 SU
 1 2 3 4 5

- 1- **BKH:** Basket Type
- 2- **5:** Catheter Diameter
- 3- **20:** Basket Opening
- 4- **12:** Catheter Length
- 5- **SU:** Handle Feature

Basket Type	BKH BKO BKS	Hexagonal Basket Catheter Oval Basket Catheter Spiral Basket Catheter
Catheter Diameter	5 7	5 F (1,8 mm) 7 F (2,3 mm)
Basket Opening	20 25 30 35	20 mm 25 mm 30 mm 35 mm
Catheter Length	12 18 21 23 25 31	120 cm 180 cm 215 cm 230 cm 250 cm 310 cm
Handle Feature	SU SY RD RY	Fixed End Fixed End Flushing Channel Rotatable Rotatable Flushing Channel